

---

# synedra Web

Feature overview  
Version 3.1.0.x "Artemis"

© 2012

synedra information technologies

[www.synedra.com](http://www.synedra.com)

## Table of Contents

---

1. Feature overview .....	1
1.1. General features .....	1
1.2. Data security .....	1
1.3. Document shares .....	1
1.4. Document retrieval for external users .....	2

## 1. Feature overview

### 1.1. General features

- Web portal for the distribution of multimedia documents (reports, image data, video sequences,...) to doctors and institutions outside of the hospital
- Customizable design options (the customer's logo can be inserted)

### 1.2. Data security

- Shares can be made available to all users of an organization unit or to just one external user
- Login required for internal and external users
- Configurable validity period for passwords
- Automatic expiration date for each shared document
- Cancellation of shares possible
- Security protocols
- Deactivation of users
- Auditing of file accesses
- Server-sided data encryption, HTTPS protocol
- The application can be run on physically separated servers

### 1.3. Document shares

- Automatic shares via HL7 communication, e.g. for documents from the HIS/RIS
- Manual sharing by internal users in synedra View and synedra Control

- Additional text messages for external users possible
- Configurable information in e-mails
- Each document share can be made available individually to exactly one user account

## 1.4. Document retrieval for external users

- External recipients receive an e-mail notification as soon as new documents have been shared
- Contact form for external users can be displayed and masked, e.g. for requesting reports
- User-specific login enables an easy and secure access via Internet
- Preview function for image and video data
- Download function, optional integration with synedra View